


Flexible Iris Retractors




- Single-use or reusable
- Strong yet flexible polypropylene hook for mechanical iris dilation
- Vertical orientation for easy handling
- Insertions through 0.5 mm self-sealing incision at the limbus
- Supplied 5 retractors per set

Round Iris Retractor



7001	Iris Retractors, 1 set of 5 retractors
7014	Iris Retractors, 1 set of 4 retractors
7002	Iris Retractors, 6 sets of 5 retractors


Flat Iris Retractor



7004	Flat Iris Retractors, 1 set of 5 retractors
7013	Flat Iris Retractors, 6 sets of 5 retractors

- Flat profile polypropylene fiber
- Less Iris trauma due to broader contact to pupil
- No twisting due to flat profile

Reusable Iris Retractor



7003	Reusable Iris Retractors, 1 set of 6 retractors
------	---

Autoclavable container

Optalmed

1249 Quarry Lane, Suite 140
Pleasanton, CA 94566 USA

p. +1.925.462.6116
f. +1.213.270.9334

info@optalmed.com
www.optalmed.com

